



CircuFlow™ 5150

Lymphedema Compression Therapy

Product Overview

The **CircuFlow™ 5150** is a lymphedema compression therapy pump used to help relieve the symptoms associated with chronic lymphatic and venous disorders in the extremities. When the device is powered on, the attached garments use sequential air compression to move lymphatic fluid throughout the body, mimicking the natural flow. This movement can reduce swelling.

Features and Benefits:

- ✓ **Lightweight:** 5.5 pounds, easy to move as needed.
- ✓ **Comfortable:** Quiet inflation and deflation, relaxing, effective therapy.
- ✓ **Safe:** Automatic shut off and emergency stop feature.

Device Specifications:

Weight:	5.5 lbs.
Inflation Cycle:	60 seconds
Vacuum Modes:	Continuous or intermittent
Rest Time:	12 seconds
Dimensions:	4.5" (H) x 9.5" (D) x 10" (W)
Cycle Time:	15 seconds / chamber
End Therapy Vacuum Time:	180 seconds
Treatment Time:	0.5 - 3 hours
Pressure Range:	20 - 80 mmHg
Electrical:	120VAC, 1.0AMP, SLO-BLO
Fuse Rated:	250VAC, 1.0AMP, SLO-BLO
Operation Mode:	Intermittent Gradient Sequential



Item Code: P-5150
LB34.0002 Rev D 20160115