




## extriCARE® 2400

### Negative Pressure Wound Therapy System

## Product Overview

The **extriCARE® 2400** pump is a lightweight, portable, battery-powered negative pressure wound therapy pump. This device, used in conjunction with our anatomically fitted line of wound dressings, may promote wound healing through the drainage and removal of wound exudates, infectious material, and tissue debris from the wound bed using continuous and or intermittent negative pressure.

## Features and Benefits:

- 
**Effective:** Powerful negative pressure, may promote faster wound healing
- 
**Reliable:** Rechargeable battery, 30 hours of continuous battery life with a full charge
- 
**Safe:** Airtight, disposable canisters that come in two sizes, gelling agent solidifies wound exudates

## Device Specifications:

<b>Weight:</b>	8.6 oz.
<b>Battery Type:</b>	Lithium (rechargeable)
<b>Vacuum Modes:</b>	Continuous or intermittent
<b>Charging Time:</b>	< 3.5 hours
<b>Dimensions:</b>	Length: 3.35" (8.5 cm) Height: 5.67" (14.4 cm) Width: 1.46" (3.7 cm)
<b>Input:</b>	100-240 V 50/60 Hz 0.4 A



Item Code: 2400  
LB30.0002 Rev E 20151006